

Product datasheet for **TP303916M**

GPBP1L1 (NM_021639) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human GC-rich promoter binding protein 1-like 1 (GPBP1L1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203916 protein sequence Red =Cloning site Green =Tags(s)
	MAQHDFVPAWLNFTSPQSAKSPTATFEKHGEHLPRGEGRFGVSRRRHNSSDGFFNNGPLRTAGDSWH QPS LFRHDSVDSGVSKGAYAGITGNPSGWHSSSRGHDGMSQRSGGGTGNHRHWNGSFHSRKGAFCQEKPP MEI REEKKEDKVEKLQFEEEDFPSLNPEAGKQHQPCTPIGTPSGVWENPPSAKQPSKMLVIKKVSKEDPAAAF SAAFTSPGSHHANGNKLSSVPSVYKNLVPKPVPPPSKPNAWKANRMEHKSGSLSSSRSAFTSPISVTK PVLASGAALSSPKESPSSTTPPIEISSRLTKLTRTTDRKSEFLKTLKDDRNGDFSENDRDCDKLEDLE DNSTPEPKENGEEGCHQNGALPVEEGEVLSSHLEAEHRLKAMGWQEYPENDENCLPLTEDELKEFH M KTEQLRRNGFGKNGFLQSRSSSLFSPWRSTCKAEFEDSDTETSSSETSDDDAWK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	52.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_067652</u>
Locus ID:	60313
UniProt ID:	<u>Q9HC44</u>
RefSeq Size:	3754
Cytogenetics:	1p34.1
RefSeq ORF:	1422
Synonyms:	SP192
Summary:	Possible transcription factor.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified GPBP1L1 protein (Cat# [TP303916]). The protein was produced from HEK293T cells transfected with GPBP1L1 cDNA clone (Cat# [RC203916]) using MegaTran 2.0 (Cat# [TT210002]).