

## Product datasheet for **TP303933**

### C12orf23 (TMEM263) (NM\_152261) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 12 open reading frame 23 (C12orf23), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203933 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MNQTDKNQQEIPSYLNDEPPEGSMKDHPPQQPGMLSRVTGGIFSVTKGAVGATIGGVAVWIGGKSLEVTKT AVTTVPSMGIGLVKGGVSAVAGGVTAVGSAVWNKVPLTGKKKDKSD  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	11.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_689474</a>
Locus ID:	90488
UniProt ID:	<a href="#">Q8WUH6</a>
RefSeq Size:	3386
Cytogenetics:	12q23.3



[View online »](#)

RefSeq ORF: 348  
Synonyms: C12orf23  
Protein Families: Transmembrane

### Product images:



Coomassie blue staining of purified TMEM263 protein (Cat# TP303933). The protein was produced from HEK293T cells transfected with TMEM263 cDNA clone (Cat# [RC203933]) using MegaTran 2.0 (Cat# [TT210002]).