

Product datasheet for **TP303941**

Clusterin (CLU) (NM_203339) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human clusterin (CLU), transcript variant 2, 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC203941 protein sequence
Red=Cloning site **Green**=Tags(s)

MMKTLLLFVGLLLTWESGQVLGDQTVSDNELQEMSNQGSKYVNKEIQNAVNGVKQIKTLIEKTNEERKTL
LSNLEEAKKKKEDALNETRESETKLKELPGVCNETMMALWEECKPCLKQTCMKFYARVCRSGSGLVGRQL
EEFLNQSSPFYFWMNGDRIDSLENDRQQTHMLDVMQDHFSSRASSIIDELFQDRFFFTREPQDTHYLPFS
LPHRRPHFFFPKSRIVRSLMPFSPYEPLNFHAMFQPFLEMIHEAQQAMDIIHFHSPAFQHPPTFEIREGDD
DRTVCREIRHNSTGCLRMKDQCDKCREILSVDCSTNNPSQAKLRREDESLQVAERLTRKYNELLKSYQW
KMLNTSSLLEQLNEQFNWWSRLANLTQGEDQYYLRVTTVASHTSDSDVPSGVTEVWVKLFDSDPITVTP
VEVSRKNPKFMETVAEKALQEYRKKHREE

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 50 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


RefSeq: [NP_976084](#)



[View online »](#)

Locus ID:	1191
UniProt ID:	P10909
RefSeq Size:	3012
Cytogenetics:	8p21.1
RefSeq ORF:	1347
Synonyms:	AAG4; APOJ; CLI; KUB1; MGC24903; SGP-2; SGP2; SP-40; TRPM-2; TRPM2
Summary:	The protein encoded by this gene is a secreted chaperone that can under some stress conditions also be found in the cell cytosol. It has been suggested to be involved in several basic biological events such as cell death, tumor progression, and neurodegenerative disorders. Alternate splicing results in both coding and non-coding variants.[provided by RefSeq, May 2011]
Protein Families:	Druggable Genome, Secreted Protein

Product images:



Coomassie blue staining of purified CLU protein (Cat# TP303941). The protein was produced from HEK293T cells transfected with CLU cDNA clone (Cat# [RC203941]) using MegaTran 2.0 (Cat# [TT210002]).