

Product datasheet for **TP303950M**

PHYHIPL (NM_032439) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phytanoyl-CoA 2-hydroxylase interacting protein-like (PHYHIPL), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC203950 protein sequence Red =Cloning site Green =Tags(s)
	MEVPRLDHALNSPTSPCEEVIKNSLEAIQLCDRDGNKSQDSGIAEMEELPVPHPNIKISNITCDSFKISW EMDSKSKDRITHYFIDLNNKENKNSNKFHKHDVPTKLVAKAVPLPMTVRGHWFLSPRTEYTVAVQTASKQ VDGDYVWSEWSEIIEFCTADYSKVHLTQLLEKAEVIAGRMLKFSVFYRNQHKKEYFDYVREHHGNAMQPSV KDNSGSHGSPISGKLEGIFFCSTEFNTGKPPQDSPYGRYRFEIAAEKLFNPNTNLYFGDFCYMYTAYHY VILVIAPVGSFGDEFCKQRLPQLNSKDNKFLTCTEEDGVLVYHHAQDVILEVIYTDVPDLSLGTVAEITG HQLMSLSTANAKKDPSCKTCNISVGR
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	42.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_115815</u>



[View online »](#)

Locus ID:	84457
UniProt ID:	<u>Q96FC7</u>
RefSeq Size:	3581
Cytogenetics:	10q21.1
RefSeq ORF:	1128
Summary:	May play a role in the development of the central system.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified PHYHIPL protein (Cat# [TP303950]). The protein was produced from HEK293T cells transfected with PHYHIPL cDNA clone (Cat# [RC203950]) using MegaTran 2.0 (Cat# [TT210002]).