

Product datasheet for **TP304023M**

TRAX (TSNAX) (NM_005999) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human translin-associated factor X (TSNAX), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204023 protein sequence Red =Cloning site Green =Tags(s)
	<p>MSNKEGSGGFRKRKHDNFPHNQRREGKDVNSSSPVMLAFKSFQQELDARHDKYERLVKLSRDITVESKRT IFLLHRITSAPDMEDILTESEIKLDGVRQKIFQVAQELSGEDMHQFHRAITGLQEYVEAVSFQHFQIKTR SLISMDEINKQLIFTTEDNGKENKTPSSDAQDKQFGTWRLRVTPVDYLLGVADLTGELMRMCINSVGNNGD IDTPFEVSQFLRQVYDGFSGFIGNTGPYEVSKKLYTLKQSLAKVENACYALKVVRGSEIPKHMLADVFSVKT EMIDQEEGIS</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	32.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005990
Locus ID:	7257



[View online »](#)

UniProt ID: [Q99598](#), [A0A024R3V8](#)

RefSeq Size: 2667


Cytogenetics: 1q42.2

RefSeq ORF: 870

Synonyms: C3PO; TRAX

Summary: This gene encodes a protein which specifically interacts with translin, a DNA-binding protein that binds consensus sequences at breakpoint junctions of chromosomal translocations. The encoded protein contains bipartite nuclear targeting sequences that may provide nuclear transport for translin, which lacks any nuclear targeting motifs. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified TSNAX protein (Cat# [TP304023]). The protein was produced from HEK293T cells transfected with TSNAX cDNA clone (Cat# [RC204023]) using MegaTran 2.0 (Cat# [TT210002]).