

Product datasheet for **TP304049L**

Transaldolase 1 (TALDO1) (NM_006755) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human transaldolase 1 (TALDO1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204049 protein sequence Red =Cloning site Green =Tags(s)

MSSSPVKRQRMESALDQLKQFTTVVADTGDFHAIDEYKPQDATTNPSLILAAAQMPAYQELVEEAIAYGR
KLGGSQEDQIKNAIDKLFVLFGAIELKKIPGRVSTEV DARLSFDKDAMVARARRLIELYKEAGISKDRIL
IKLSSTWEGIQAGKELEEQHGIHCNMTLLFSFAQAVACAEAGVTLISPFVGRILDWHVANTDKKSYEPL
DPGVKSVTKIYNYYKKFSYKTIVMGASFRNTGEIKALAGCDFLTISP KLLGELLQDNAKLVVLSAKAAQ
ASDLEKIHLDEKSFRLHNE DQMAVEKLS DGIRKFAADAVKLERMLTERM FNAENGK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	37.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006746
Locus ID:	6888



[View online »](#)

UniProt ID: [P37837](#), [A0A140VK56](#)

RefSeq Size: 1319

Cytogenetics: 11p15.5

RefSeq ORF: 1011


Synonyms: TAL; TAL-H; TALDOR; TALH

Summary: Transaldolase 1 is a key enzyme of the nonoxidative pentose phosphate pathway providing ribose-5-phosphate for nucleic acid synthesis and NADPH for lipid biosynthesis. This pathway can also maintain glutathione at a reduced state and thus protect sulfhydryl groups and cellular integrity from oxygen radicals. The functional gene of transaldolase 1 is located on chromosome 11 and a pseudogene is identified on chromosome 1 but there are conflicting map locations. The second and third exon of this gene were developed by insertion of a retrotransposable element. This gene is thought to be involved in multiple sclerosis. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Pentose phosphate pathway

Product images:



Coomassie blue staining of purified TALDO1 protein (Cat# [TP304049]). The protein was produced from HEK293T cells transfected with TALDO1 cDNA clone (Cat# [RC204049]) using MegaTran 2.0 (Cat# [TT210002]).