

Product datasheet for **TP304072M**

alpha Tubulin (TUBA3C) (NM_079836) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human tubulin, alpha 3c (TUBA3C), transcript variant 2, 100 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC204072 protein sequence
Red=Cloning site Green=Tags(s)

MRECISIHVGQAGVQIGNACWELYCLEHGIQPDGQMPSDKTIGGGDDSFNTFFSETGAGKHVPRAVFVDL
EPTVVDEVRTGTYRQLFHPEQLITGKEDAANNYARGHYTIGKEIVDLVLDRIKRLADLCTGLQGFLIFHS
FGGGTGS GFASLLMERLSVDYGKSKLEFAIYPAPQVSTAVVEPYNSILTHTTLEHSDCAFMDVNEAIY
DICRRNLDIERPTYTNLNLIGIQIVSSITASLRFDGALNVDLTEFQTNLVPYPRIFPLATYAPVISA EK
AYHEQLSVAEITNACFEPANQMVKCDPRHGKYMCCMLYRGDWPVKDVNAAIATIKTKRTIQFVDWCPTG
FKVGINYQPPTVVPGGDLAKVQRAVCMLSNTTIAEAWARLDAALEKDYEEVGVDSVEAEAEEGEEY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 49.8 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Bioactivity: ELISA capture for autoantibodies (PMID: [29427009](https://pubmed.ncbi.nlm.nih.gov/29427009/))

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_524575](https://www.ncbi.nlm.nih.gov/nuccore/NP_524575)



[View online »](#)


Locus ID:	7278
UniProt ID:	Q13748
RefSeq Size:	1408
Cytogenetics:	13q12.11
RefSeq ORF:	1254
Synonyms:	bA408E5.3; TUBA2; tubulin, alpha 2; tubulin, alpha 3c

Summary: Microtubules of the eukaryotic cytoskeleton perform essential and diverse functions and are composed of a heterodimer of alpha and beta tubulin. The genes encoding these microtubule constituents are part of the tubulin superfamily, which is composed of six distinct families. Genes from the alpha, beta and gamma tubulin families are found in all eukaryotes. The alpha and beta tubulins represent the major components of microtubules, while gamma tubulin plays a critical role in the nucleation of microtubule assembly. There are multiple alpha and beta tubulin genes and they are highly conserved among and between species. This gene is an alpha tubulin gene that encodes a protein 99% identical to the mouse testis-specific Tuba3 and Tuba7 gene products. This gene is located in the 13q11 region, which is associated with the genetic diseases Clouston hidrotic ectodermal dysplasia and Kabuki syndrome. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Stem cell - Pluripotency

Protein Pathways: Gap junction, Pathogenic Escherichia coli infection

Product images:



Coomassie blue staining of purified TUBA3C protein (Cat# [TP304072]). The protein was produced from HEK293T cells transfected with TUBA3C cDNA clone (Cat# [RC204072]) using MegaTran 2.0 (Cat# [TT210002]).