

Product datasheet for TP304079

AIPL1 (NM_014336) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human aryl hydrocarbon receptor interacting protein-like 1 (AIPL1), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204079 protein sequence Red=Cloning site Green=Tags(s)

MDAALLNVEGVKKTILHGGTGELPNFITGSRVIFHFRTMKCDEERTVIDDSRQVGQPMHIIIGNMFKLE
VWEILLTSMRVHEVAEFWCDTIHTGVYPILSRSLRQMAQGKDPTTEWHVHTCGLANMFAYHTLGVEDLDEL
QKEPQPLVFIELLQVDAPSDYQRETWNLSNHEKMKAVPVLHGEGNRLFGLGRYEEASSKYQEAIICLRN
LQTKKEPWVQWLKLEKMINTLILNYCQCLLKKEEYEVLEHTSDILRHHPGIVKAYYVRARAHAEVWNE
AEAKADLQKVLELEPSMQKAVRRELRLLENRMAEKQEEERLCRNMLSQGATQPPAEPPTPEPPAQSTEP
PAEPPTAPSAELSAGPPAEPATEPPPPSPGHSLQH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	43.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_055151</u>




[View online »](#)

Locus ID: 23746
UniProt ID: [Q9NZN9](#), [F1T0B6](#)
RefSeq Size: 2990
Cytogenetics: 17p13.2
RefSeq ORF: 1152
Synonyms: AIPL2; LCA4

Summary: Leber congenital amaurosis (LCA) is the most severe inherited retinopathy with the earliest age of onset and accounts for at least 5% of all inherited retinal diseases. Affected individuals are diagnosed at birth or in the first few months of life with nystagmus, severely impaired vision or blindness and an abnormal or flat electroretinogram. The photoreceptor/pineal-expressed gene, AIPL1, encoding aryl-hydrocarbon interacting protein-like 1, is located within the LCA4 candidate region. The encoded protein contains three tetratricopeptide motifs, consistent with chaperone or nuclear transport activity. Mutations in this gene may cause approximately 20% of recessive LCA. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified AIPL1 protein (Cat# TP304079). The protein was produced from HEK293T cells transfected with AIPL1 cDNA clone (Cat# [RC204079]) using MegaTran 2.0 (Cat# [TT210002]).