

Product datasheet for **TP304095**

SURF2 (NM_017503) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human surfeit 2 (SURF2), 20 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC204095 protein sequence
Red=Cloning site **Green**=Tags(s)

MSELPGDVRAFLREHPSRLQTDARKVRCILTGHELPCRLPELQVYTRGKKYQRLVRASPAFDYAEFEPH
IVPSTKNPHQLFCKLTLRHINKCPEHVLRLHTQGRRYQALCKYECCQKQGVYVPACLVHRRRRREDQMD
GDGPRPREAFWEPTSSDEGGAASDDSM TDLYPPEL FTRKDLGSTEDGDGTDDFLTDKEDEKAKPPREKAT
DEGRRETTVYRGLVQKRGKKQLGSLKKKFKSHHRKPKSFSSCKQPG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 29.5 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP_059973

Locus ID: 6835

UniProt ID: Q15527



[View online »](#)

RefSeq Size: 851


Cytogenetics: 9q34.2

RefSeq ORF: 768

Synonyms: SURF-2

Summary: This gene shares a bidirectional promoter with surfeit 1 (SURF1; GeneID: 6834), which is located on the opposite strand. It encodes a conserved protein that is expressed in a variety of tissues. [provided by RefSeq, Jul 2013]

Product images:



Coomassie blue staining of purified SURF2 protein (Cat# TP304095). The protein was produced from HEK293T cells transfected with SURF2 cDNA clone (Cat# [RC204095]) using MegaTran 2.0 (Cat# [TT210002]).