

Product datasheet for **TP304119**

PTCD3 (NM_017952) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human Pentatricopeptide repeat domain 3 (PTCD3), 20 µg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC204119 representing NM_017952
Red=Cloning site **Green**=Tags(s)

MAVVS AVRWLGLRSRLGQPLTGRRAGLCEQARSCR FYSGSATLSKVEGTDVTGIEEVVIPKKKTW DQKVAV
LQALASTVNRD TTAVPYVFQDDPYLMPASSLESRSFLLAKKSGEN VAKFIINSYPKYFQKDIAEPHIPCL
MPEYFEPQIKDISEAALKERIELRKVKASVDMFDQLLQAGTTVSLETTNSLLDLLCCYGDQEPSTDYHFQ
QTGQSEALEEENDETSRRKAGHQFGVTWRAKNNAERIFSLMPEKNEHSYCTMIRGMV KHRAEQALNLYT
ELLNNR LHADVYTFNALIEATVCAINEKFEEKWSKILELLRHMVAQKV KPNLQTFNTILKCLRRFHFVFAR
SPALQVLREMKAIGIEPSLATYHHIIRLFDQPGDPLKRSSFIIYDIMNELMGKRFSPKDPDDDKFFQSAM
SICSSLRDLELAYQVHG LLKTGDNWKFIGPDQHRNFYSKFFDLICLMEQIDVTLKWYEDLIPSAYFPHS
QTMIHLLQALDVANRLEVIPKIWKDSKEYGHTFRSDLREEILMLMARDKHPELQVAFADCAADIKSAYE
SQPIRQTAQDWPATSLNCIAILFRAGRTQEAWKMLGLFRKHNKIPRSELLNELMDSAKVSN SPSQAIEV
VELASAFSLPICEGLTQRVMSDFAINQEQKEALS NLTALTS DSDTSSSDSDSDTSE GK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag: C-Myc/DDK
Predicted MW: 78.4 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage: Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_060422
Locus ID:	55037
UniProt ID:	Q96EY7
RefSeq Size:	6708
Cytogenetics:	2p11.2
RefSeq ORF:	2067
Synonyms:	MRP-S39
Summary:	Mitochondrial RNA-binding protein that has a role in mitochondrial translation. [UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified PTCD3 protein (Cat# TP304119). The protein was produced from HEK293T cells transfected with PTCD3 cDNA clone (Cat# [RC204119]) using MegaTran 2.0 (Cat# [TT210002]).