

Product datasheet for TP304156

Superoxide Dismutase 3 (SOD3) (NM_003102) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human superoxide dismutase 3, extracellular (SOD3), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204156 protein sequence Red =Cloning site Green =Tags(s)
	MLALLCSCLLLAAGASDAWTGEDSAEPNSDSA EWIRDMYAKVTEIWQEV MQRRDDDGLHAACQVQPS AT LDAAQPRVTGVVLFRLAPRAKLDAFFALEGFTEPNSSSRAIHVHQFGDLSQGCESTGPHYNPLAVPHP QHPGDFGNFAVRDGLWRYRAGLAASLAGPHSIVGRAV VVHAGEDDLGRGGNQASVENGNAGRRLACC VV GVCGPGLWERQAREHSERKKRRRESECKAA
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	24.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003093
Locus ID:	6649



[View online »](#)

UniProt ID: [P08294](#)

RefSeq Size: 1546

Cytogenetics: 4p15.2


RefSeq ORF: 720

Synonyms: EC-SOD

Summary: This gene encodes a member of the superoxide dismutase (SOD) protein family. SODs are antioxidant enzymes that catalyze the conversion of superoxide radicals into hydrogen peroxide and oxygen, which may protect the brain, lungs, and other tissues from oxidative stress. Proteolytic processing of the encoded protein results in the formation of two distinct homotetramers that differ in their ability to interact with the extracellular matrix (ECM). Homotetramers consisting of the intact protein, or type C subunit, exhibit high affinity for heparin and are anchored to the ECM. Homotetramers consisting of a proteolytically cleaved form of the protein, or type A subunit, exhibit low affinity for heparin and do not interact with the ECM. A mutation in this gene may be associated with increased heart disease risk. [provided by RefSeq, Oct 2015]

Protein Families: Druggable Genome, Secreted Protein

Product images:



Coomassie blue staining of purified SOD3 protein (Cat# TP304156). The protein was produced from HEK293T cells transfected with SOD3 cDNA clone (Cat# [RC204156]) using MegaTran 2.0 (Cat# [TT210002]).