

Product datasheet for **TP304219M**

CCDC97 (NM_052848) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 97 (CCDC97), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204219 protein sequence Red =Cloning site Green =Tags(s)
	<p>MEAVATATAAKEPDKGCIEPGPGHWGELSRTVPVSKPQDKVEAAEATPVALDSDTSGAENAAVSAMLHAV AASRLPVCSQQQGEPDLTEHEKVAILAQLYHEKPLVFLERFRTGLREEHLACFGHVGRDHRADFYCAEVA RQGTARPRTLRLRNRRYAALRELIQGGYFSDEQMRFRAPLLYEQYIGQYLQEELSARTPTHQPPKP GSPGRPACPLSNLLLSYEERELQQRLLQQEEEEACLEEEEEEDSDEEDQRSGKDSEAWVPDSEERLI LREEFTSRMHQRFLDGKDGDFDYSTVDDNPDFDNLDIRARDEEERYFDEEEPEDAPSPELDGD</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	38.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_443080</u>
Locus ID:	90324



[View online »](#)

UniProt ID: [Q96F63](#)
RefSeq Size: 3344
Cytogenetics: 19q13.2
RefSeq ORF: 1029

Product images:



Coomassie blue staining of purified CCDC97 protein (Cat# [TP304219]). The protein was produced from HEK293T cells transfected with CCDC97 cDNA clone (Cat# [RC204219]) using MegaTran 2.0 (Cat# [TT210002]).