

Product datasheet for **TP304244M**

SEPTIN6 (NM_145802) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human septin 6 (SEPT6), transcript variant V, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204244 protein sequence Red =Cloning site Green =Tags(s)
	<p>MAATDIARQVGEGRCTVPLAGHVGFDLSPDQLVNKSVSQGFCFNILCVGETGLGKSTLMDTLFNTKFEGE PATHTQPGVQLQSNTYDLQESNVRLKLTIVSTVGFQDQINKEDSYKPIVEFIDAQFEAYLQEELKIRRVL HTYHDSRIHVCLYFIAPTGHSLKSLDLVTMKKLDSKVNIIIAKADAISKSELTKFKIKITSELVSNV QIYQFPTDDESVAEINGTMNAHLPFAVIGSTEELKIGNKMMRARQYPWGTQVQVENEAHCDFVKLREMLIR VNMEDLREQTHTRHYELRRCLEEMGFKDTPDPSKPFSLQETYEAKRNEFLGELQKKEEEMRQMFVQRV KEKEAELKEAKELHEKFDRLKHLHQDEKKKLEDKKKSLDDEVNAFKQRKTA AELLQSQGSQAGGSQTLK RDKEKKNFF</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	49 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_665801</u>




[View online »](#)

Locus ID: 23157
UniProt ID: [Q14141](#), [Q548C9](#)
RefSeq Size: 2564
Cytogenetics: Xq24
RefSeq ORF: 1287
Synonyms: SEP2; SEPT2; SEPT6

Summary: This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis. One version of pediatric acute myeloid leukemia is the result of a reciprocal translocation between chromosomes 11 and X, with the breakpoint associated with the genes encoding the mixed-lineage leukemia and septin 2 proteins. This gene encodes four transcript variants encoding three distinct isoforms. An additional transcript variant has been identified, but its biological validity has not been determined. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified SEPTIN6 protein (Cat# [TP304244]). The protein was produced from HEK293T cells transfected with SEPTIN6 cDNA clone (Cat# [RC204244]) using MegaTran 2.0 (Cat# [TT210002]).