

Product datasheet for **TP304248M**

NKIRAS1 (NM_020345) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human NFkB inhibitor interacting Ras-like 1 (NKIRAS1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204248 protein sequence Red =Cloning site Green =Tags(s)
	MGKGCKVWVCGLLSVGKTAILEQLLYGNHTIGMEDCETMEDVYMASVETDRGVKEQLHLHYDTRGLQEGVE LPKHYSFADGFVLVYSVNNLESFQRVELLKKEIDKFKDKKEVAIWLGNKIDLSEQRQVDAEVAQQWAK SEKVRLEWTVTDRKTLIEPFTLLASKLSQPQSKSSFPLPGRKNKGNSNSEN
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	21.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_065078</u>
Locus ID:	28512
UniProt ID:	<u>Q9NYS0</u>
RefSeq Size:	2013



[View online »](#)

Cytogenetics:	3p24.2
RefSeq ORF:	576
Synonyms:	kappaB-Ras1; KBRAS1
Summary:	<p>Atypical Ras-like protein that acts as a potent regulator of NF-kappa-B activity by preventing the degradation of NF-kappa-B inhibitor beta (NFKBIB) by most signals, explaining why NFKBIB is more resistant to degradation. May act by blocking phosphorylation of NFKBIB and mediating cytoplasmic retention of p65/RELA NF-kappa-B subunit. It is unclear whether it acts as a GTPase. Both GTP- and GDP-bound forms block phosphorylation of NFKBIB.</p> <p>[UniProtKB/Swiss-Prot Function]</p>

Product images:



Coomassie blue staining of purified NKIRAS1 protein (Cat# [TP304248]). The protein was produced from HEK293T cells transfected with NKIRAS1 cDNA clone (Cat# [RC204248]) using MegaTran 2.0 (Cat# [TT210002]).