

## Product datasheet for **TP304263**

### **C7orf30 (MALSU1) (NM\_138446) Human Recombinant Protein**

#### **Product data:**

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 7 open reading frame 30 (C7orf30), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204263 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MGPGRVARLLAPLMWRRVSSVAGSAVGAEPGLRLLAVQRLPVGAAFCRACQTPNFVRGLHSEPGLEER</p> <p>AEGTVNEGRPESDAADHTGPKFDIDMMVSLRQENARDICVIQVPPEMRYTDYFVIVSGTSTRHLHAMAF</p> <p>YVVKMYKHLKCKRDPHVKIEGKDTDDWLCVDFGSMVIHLMLPETREIYELEKLWTLRSYDDQLAQIAPET</p> <p>VPEDFILGIEDDTSSVTPVELKCE</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	26 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_612455</a></u>
Locus ID:	115416



[View online »](#)

UniProt ID: [Q96EH3](#)

RefSeq Size: 774


Cytogenetics: 7p15.3

RefSeq ORF: 702

Synonyms: C7orf30; mtRsfA

**Summary:** Required for normal mitochondrial ribosome function and mitochondrial translation (PubMed:22238375, PubMed:23171548). May play a role in ribosome biogenesis by preventing premature association of the 28S and 39S ribosomal subunits (Probable). Interacts with mitochondrial ribosomal protein L14 (MRPL14), probably blocking formation of intersubunit bridge B8, preventing association of the 28S and 39S ribosomal subunits (Probable). Addition to isolated mitochondrial ribosomal subunits partially inhibits translation, probably by interfering with the association of the 28S and 39S ribosomal subunits and the formation of functional ribosomes (Probable). May also participate in the assembly and/or regulation of the stability of the large subunit of the mitochondrial ribosome (PubMed:22238376, PubMed:23171548). May function as a ribosomal silencing factor (Probable).[UniProtKB/Swiss-Prot Function]

## Product images:



Coomassie blue staining of purified MALSU1 protein (Cat# TP304263). The protein was produced from HEK293T cells transfected with MALSU1 cDNA clone (Cat# [RC204263]) using MegaTran 2.0 (Cat# [TT210002]).