

Product datasheet for **TP304267**

GMP Synthase (GMPS) (NM_003875) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human guanine monophosphate synthetase (GMPS), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204267 protein sequence Red =Cloning site Green =Tags(s)

MALCNGDSKLENAGGDLKDGHHHYEGAVILDAGAQQYKVIDRRVRELFVQSEIFPLETPAFAIKEQGFR
AIIISGGPNSVYAEDAPWFDPAIFTIGKPVLGICYGMQMMNKVFGGTVHKKSVDREDFVFNISVDNTCSLF
RGLQKEEVLLTHGDSVDKVADGFKVARSIGNIVAGIANESKKLYGAQFHPEVGLTENGKVLKNFLYDI
AGCSGFTTVQNRELECIKERVGTSKVLVLLSGGVDSTVCTALLNRALNQEVIQVIHIDNGFMRKRES
QSVEEALKKLGIVKVINAAHSFYNGTTTLPISDEDRTPRKRISKTLNMTTSPEEKRKIIIGDTFVKIANE
VIGEMNLKPEEVFLAQGTLRPDLIESASLVASGKAEIKTHHNDTELIRKLREEGKVIEPLKDFHKDEVR
ILGRELGLPEELVSRHPFPGPLAIRVICAEEPYICKDFPETNNILKIVADFSASVKKPHTLLQRVKACT
TEEDQEKLMTSLHSLNAFLLPKTVGVQGDRCRSYSYVCGISSKDEPDWESLIFLARLIPRMCHNVNRV
VYIFGPPVKEPPTDVTPTFLTTGVLSTLRQADFEAHNILRESGYAGKISQMPVILTPHFDRDPLQKQPS
CQRSVIRTFITSDFMGTGIPATPGNEIPVEVLKVMTEIKKIPGISRIMYDLTSKPPGTTEWE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	76.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_003866](#)

Locus ID: 8833

UniProt ID: [P49915](#), [A0A140VJK6](#)

RefSeq Size: 2457

Cytogenetics: 3q25.31

RefSeq ORF: 2079


Synonyms: GATD7

Summary: In the de novo synthesis of purine nucleotides, IMP is the branch point metabolite at which point the pathway diverges to the synthesis of either guanine or adenine nucleotides. In the guanine nucleotide pathway, there are 2 enzymes involved in converting IMP to GMP, namely IMP dehydrogenase (IMPD1), which catalyzes the oxidation of IMP to XMP, and GMP synthetase, which catalyzes the amination of XMP to GMP. [provided by RefSeq, Jul 2008]

Protein Families: Stem cell - Pluripotency

Protein Pathways: Drug metabolism - other enzymes, Metabolic pathways, Purine metabolism

Product images:



Coomassie blue staining of purified GMPS protein (Cat# TP304267). The protein was produced from HEK293T cells transfected with GMPS cDNA clone (Cat# [RC204267]) using MegaTran 2.0 (Cat# [TT210002]).