

## Product datasheet for TP304325

### beta glucuronidase (GUSB) (NM\_000181) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glucuronidase, beta (GUSB), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204325 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MARGSAVAWAALGPLLWGCALGLQGGMLYPQESPSRECKELDGLWSFRADFSNRRRGFEEQWYRRPL  
WE  
SGPTVDMPVPSSFNDISQDWRLRHVFGWVWYEREVILPERWTQDLRTRVVLRIGSAHSYAIWVNGVDTL  
EHEGGYLPFEADISNLVQVGPLPSRLRITIAINNTLPTTLPPGTIQYLTDTSKYPKGYFVQNTYFDFFN  
YAGLQRSVLLYTTPTYIDDITVTSVEQDSGLVNYQISVKGSNLFKLEVRLLDAENKVVANGTGTQGQL  
KVPGVSLWWPYLMHERPAYLYSLEVQLTAQTSGLPVSDFYTLVPGIRTVAVTKSQFLINGKPFYFHGVNK  
HEDADIRGKGFDPVLLVKDFNLLRWLGANAFRTSHYPYAEEMQMCDRYGIWIDECPGVGLALPQFFN  
N  
VSLHHHMQVMEEVRRDKNHPAVVMWSVANEPASHLESAGYYLKMVIAHTKSLDPSRPVTFVSNSNYA  
AD  
KGAPYVDVICLNSYYSWYHDYGHLELIQLQLATQFENWYKKYQKPIIQSEYGAETIAGFHQDPPLMFTEE  
YQKSLEQYHLGLDQRRKYVWGELIWNFADFMTQSPTRVLGNKKGIFTRQRQPKSAAFLLRERYWKIA  
NETRYPHSVAKSQCLENSLFT

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**


Tag:	C-Myc/DDK
Predicted MW:	72.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_000172</a>
<b>Locus ID:</b>	2990
<b>UniProt ID:</b>	<a href="#">P08236</a>
<b>RefSeq Size:</b>	2321
<b>Cytogenetics:</b>	7q11.21
<b>RefSeq ORF:</b>	1953
<b>Synonyms:</b>	BG; MPS7
<b>Summary:</b>	This gene encodes a hydrolase that degrades glycosaminoglycans, including heparan sulfate, dermatan sulfate, and chondroitin-4,6-sulfate. The enzyme forms a homotetramer that is localized to the lysosome. Mutations in this gene result in mucopolysaccharidosis type VII. Alternative splicing results in multiple transcript variants. There are many pseudogenes of this locus in the human genome.[provided by RefSeq, May 2014]
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Protein Pathways:</b>	Drug metabolism - other enzymes, Glycosaminoglycan degradation, Lysosome, Metabolic pathways, Pentose and glucuronate interconversions, Porphyrin and chlorophyll metabolism, Starch and sucrose metabolism

### Product images:



Coomassie blue staining of purified GUSB protein (Cat# TP304325). The protein was produced from HEK293T cells transfected with GUSB cDNA clone (Cat# [RC204325]) using MegaTran 2.0 (Cat# [TT210002]).