

Product datasheet for **TP304360L**

HNRPAB (HNRNPAB) (NM_031266) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human heterogeneous nuclear ribonucleoprotein A/B (HNRNPAB), transcript variant 1, 1 mg

Species: Human

Expression Host: HEK293T

Expression cDNA >RC204360 representing NM_031266

Clone or AA **Red**=Cloning site **Green**=Tags(s)

Sequence:

MSEAGEEQPMETTGTATENGHEAVPEGESPAGAGTGAAAGAGGATAAPPSGNQNGAEGDQINASKNEEDAG
KMFVGGLSWDTSKKDLKDYFTKFGVVDCTIKMDPNTGRSRGFGFILFKDAASVEKVLDQKEHRLDGRVI
DPKKAMAMKKDPVKKIFVGGNPEATEEKIREYFGEFGEIEAIELPMDPKLNKRRGFVITFKEEEPVKK
VLEKKFHTVSGSKCEIKVAQPKEVYQQQQYSGGGRGNRNRGNRSGGGGGGGGQSQSWNQGYGNYWNQGY
GYQQGYGPGYGGYDYSYPYGYGYGPGYDYSQGSTNYGKSQRRGGHQNNYKPY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 35.8 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_112556](#)

Locus ID: 3182



[View online »](#)

UniProt ID: [Q99729](#)

RefSeq Size: 1837


Cytogenetics: 5q35.3

RefSeq ORF: 996

Synonyms: ABBP1; HNRNPAB

Summary: This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are produced by RNA polymerase II and are components of the heterogeneous nuclear RNA (hnRNA) complexes. They are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene, which binds to one of the components of the multiprotein editosome complex, has two repeats of quasi-RRM (RNA recognition motif) domains that bind to RNAs. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified HNRNPAB protein (Cat# [TP304360]). The protein was produced from HEK293T cells transfected with HNRNPAB cDNA clone (Cat# [RC204360]) using MegaTran 2.0 (Cat# [TT210002]).