

Product datasheet for **TP304362L**

SHC (SHC1) (NM_003029) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human SHC (Src homology 2 domain containing) transforming protein 1 (SHC1), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204362 protein sequence Red =Cloning site Green =Tags(s) MNKLSGGGRRRTRVEGGQLGGEWTRHGSFVNKPTRGWLHPNDKVMGPGVSYLVRYMGCVEVLQSMRALD FNTRTQVTREAI SLVCEAVPGAKGATRRRKPCSRPLSSILGRSNLKFAGMPITLTVSTSSLNLMAADCKQ IIANHHMQSISFASGGDPDTAEYVAYVAKDPVNQRACHILECPEGLAQDVISTIGQAFELRFKQYLRNPP KLVTPHDRAMAGFDGSAWDEEEEEPPDHQYYNDFPGKEPPLGGVDMRLREGAAPGAARPTAPNAQTPSHL GATLPVGQPVGGDPEVRKQMPPPPPCPAGRELFDDPSYVNVQNLDKARQAVGGAGPPNPAINGSAPRDLF DMKPFEDALRVPPPQSVSMAEQLRGEPWFHKGKLSRREAEALLQLNGDFLVRESTTTPGQYVLTGLQSGQ PKHLLLVDPEGVVRTKDHRESVSHLISYHMDNHLPIISAGSELCLQPVERKL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	51.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq: [NP_003020](#)

Locus ID: 6464

UniProt ID: [P29353](#)

RefSeq Size: 3195

Cytogenetics: 1q21.3

RefSeq ORF: 1422


Synonyms: SHC; SHCA

Summary: This gene encodes three main isoforms that differ in activities and subcellular location. While all three are adapter proteins in signal transduction pathways, the longest (p66Shc) may be involved in regulating life span and the effects of reactive oxygen species. The other two isoforms, p52Shc and p46Shc, link activated receptor tyrosine kinases to the Ras pathway by recruitment of the GRB2/SOS complex. p66Shc is not involved in Ras activation. Unlike the other two isoforms, p46Shc is targeted to the mitochondrial matrix. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2011]

Protein Families: Druggable Genome

Protein Pathways: Adherens junction, Arrhythmogenic right ventricular cardiomyopathy (ARVC), Chemokine signaling pathway, Chronic myeloid leukemia, Dilated cardiomyopathy, ErbB signaling pathway, Focal adhesion, Glioma, Hypertrophic cardiomyopathy (HCM), Insulin signaling pathway, Leukocyte transendothelial migration, Natural killer cell mediated cytotoxicity, Neurotrophin signaling pathway, Pathogenic Escherichia coli infection, Regulation of actin cytoskeleton, Tight junction, Vibrio cholerae infection, Viral myocarditis

Product images:



Coomassie blue staining of purified SHC1 protein (Cat# [TP304362]). The protein was produced from HEK293T cells transfected with SHC1 cDNA clone (Cat# [RC204362]) using MegaTran 2.0 (Cat# [TT210002]).