

Product datasheet for **TP304363L**

EVI5L (NM_145245) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ecotropic viral integration site 5-like (EVI5L), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204363 protein sequence Red =Cloning site Green =Tags(s)

MASPTLSPDSSSQEALSAPTCSPTSDSENLSPELELLAKLEEQRRLLEADSKSMRSMNGSRRNSGSSLV
SSSSASSNLSHLEEDTWILWGRIANEWEWRRRKEKLLKELIRKGIPHHFRAIVWQLLCSATDMPVKNQY
SELLKMSSPCEKLIRRDIAARTYPEHEFFKGQDSLGEVLFNVMKAYSLVDREVGVCQGSFIVGLLLMQM
PEEEAFVFRMLMQEYRLRELFKPSMAELGLCIYQFEYMLQEQLPDLNTHFRSQSFHTSMYASSWFLTLF
LTTFPLPVATRVDFIMYEGLEIVFRVGLALLQVNQAELMQLDMEGMSQYFQRVIPHQFDSCPKLVKA
YQVKYNPKMKRLEKEYAAMKSKEMEEQIEIKRLRTENRLLKQRIETLEKGQVTRAQEAENYVIKRELA
VVRQQCSSAAEDLQKAQSTIRQLQEQQENPRLTEDFVSHLETELEQSRLRETETLGALREMQDKVLDMEK
RNSSLPDENNVAAQLQEELKALKVREGQAVASTRELKLQLQELSDTWQAHARGGRWKESPRKLVVDELQ
D
ELMSVRLREAQALAEGRRLRQRVELETQDHIHRNLLNRVEAERAALQEKLYLAAQNKGLQTQLSESRR
KQAEAECKSKEEVMVRLREADSMAVAEMRQRIAELEIQREEGRIQGQLNHSQYIRELKDQIEELK
AEVRLKGPFPFEDPLAFDGLSLARHLDEDSLPSSEELLGVGVGAALQDALYPLSPRDARFFRLRPA
KDSEGSSDSDADELAAPYSQGLDN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	91.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



View online »

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_660288
Locus ID:	115704
UniProt ID:	Q96CN4
RefSeq Size:	3834
Cytogenetics:	19p13.2
RefSeq ORF:	2382
Summary:	Functions as a GTPase-activating protein (GAP) with a broad specificity.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified EVI5L protein (Cat# [TP304363]). The protein was produced from HEK293T cells transfected with EVI5L cDNA clone (Cat# [RC204363]) using MegaTran 2.0 (Cat# [TT210002]).