

## Product datasheet for **TP304379M**

### ZNF385 (ZNF385A) (NM\_015481) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger protein 385A (ZNF385A), transcript variant 3, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204379 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MQPPLDLKQILPFPLEPAPTLGLFSNYSTMDPVQKAVLSHTFGGPLLTKRQVISCNICQIRFNSQSQAE  
AHYKGNRHARRVKGIEAAKTRGREPGVREPGDPAPPGSTPTNGDGVAPRPVSMENGLGPAGSPEKQPGS  
PSPSIPETGQGVTKGEGGTPAPASLPGGSKEEEEKAKRLLYCALCKVAVNSLSQLEAHNKGTKHKHILE  
ARSLGPIKAYPRLGPPTGPEPEAPAQDRTFHCEICNVKVNSEVQLKQHISSRRHRDGVAGKPNLLSRH  
KKSARGAGELAGTLTFSKELPKSLAGLLPSPLAVAAVMAAAAGSPLSLRPAPAAPLLQGPPIHPLLHPA  
PGPIRTAHGPILFSPY

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**


Tag:	C-Myc/DDK
Predicted MW:	38.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_056296</a>
Locus ID:	25946



[View online »](#)

UniProt ID:	<a href="#">Q96PM9</a>
RefSeq Size:	2444
Cytogenetics:	12q13.13
RefSeq ORF:	1098
Synonyms:	HZF; RZF; ZFP385; ZNF385
Summary:	Zinc finger proteins, such as ZNF385A, are regulatory proteins that act as transcription factors, bind single- or double-stranded RNA, or interact with other proteins (Sharma et al., 2004 [PubMed 15527981]).[supplied by OMIM, Oct 2008]
Protein Families:	Transcription Factors

### Product images:



Coomassie blue staining of purified ZNF385A protein (Cat# [TP304379]). The protein was produced from HEK293T cells transfected with ZNF385A cDNA clone (Cat# [RC204379]) using MegaTran 2.0 (Cat# [TT210002]).