

Product datasheet for TP304381

PTOP (ACD) (NM_001082486) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human adrenocortical dysplasia homolog (mouse) (ACD), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204381 protein sequence Red =Cloning site Green =Tags(s)
	<p>MPGRCQSDAAMRVNGPASRAPAGWTSGSLHTGPRAGRPRQAARGVVRGRLLLRPPPAKELPLPRKGGAWA PAGNPGPLHPLGVAVGMAGSGRLVLRPWIRELILGSETPSSPRAGQLLEVLQDAEAAVAGPSHAPDTS GATLLVSDGTHSVRCLVTREALDTSWEEKEFGFRGTEGRLLLQDCGVHVQVAEGGAPAEFYLVQVDRFS LLPTEQPRLRVPGCNQDLQKLYDCLEEHLSESTSSNAGLSLSQLLEDEMREDQEHQALVCLAESCLT LEGPCTAPPVTHWAASRCKATGEAVYTPSSMLCISENDQLLSSLGPCQRTQGPELPPDPALQDLSLT LIASPPSSPSSSGTPALPGHMSSEESGTSISLLPALSLAAPDPGQRSSSQPSAICSAPATLTPRSPHAS RTPSSPLQSCPTSLSPRSHVPSPHQALVTRPQKPSLEFKEFVGLPCKNRPPFRTGATRGAQEPCSVWEP PKRHRDGSFAFYEYEPCTSLCARVQAARLPPQLMAWALHFLMDAQPGEPTPM</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	57.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001075955
Locus ID:	65057
UniProt ID:	Q96AP0 , A0A590TQL1 , A0A5F9ZCY9
RefSeq Size:	2095
Cytogenetics:	16q22.1
RefSeq ORF:	1632
Synonyms:	PIP1; PTOP; TINT1; TPP1
Summary:	This gene encodes a protein that is involved in telomere function. This protein is one of six core proteins in the telosome/shelterin telomeric complex, which functions to maintain telomere length and to protect telomere ends. Through its interaction with other components, this protein plays a key role in the assembly and stabilization of this complex, and it mediates the access of telomerase to the telomere. Multiple transcript variants encoding different isoforms have been found for this gene. This gene, which is also referred to as TPP1, is distinct from the unrelated TPP1 gene on chromosome 11, which encodes tripeptidyl-peptidase I. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified ACD protein (Cat# TP304381). The protein was produced from HEK293T cells transfected with ACD cDNA clone (Cat# [RC204381]) using MegaTran 2.0 (Cat# [TT210002]).