

Product datasheet for **TP304392M**

PTP alpha (PTPRA) (NM_080841) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein tyrosine phosphatase, receptor type, A (PTPRA), transcript variant 3, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204392 protein sequence Red=Cloning site Green=Tags(s)

MDSWFILVLLGSLICVSANNATTVAPSVGITRLINSSTAEPVKEEAKTSNPTSSLTSLSVAPTFSPNIT
 LGPTYLTTVNSSDSNDGTTTRASTNSIGITISPNGTWLPDNQFTDARTEPWEGNSSTAATTPETFPSPDE
 TPIIAVMVALSLLVIVFIIIVLYMLRFKKYKQAGSHSNSFRLSNGRTEDVEPQSVPLLARSPSTNRKYP
 PLPVDKLEEEINRRMADDNKLFREEFNALPACPIQATCEAASKEENKEKNRYVNILPYDHSRVHLTPVEG
 VPDSYINASFINGYQEKNKFIAAQGPKEETVNDFWRMIWEQNTATIVMVTNLKERKECKCAQYWPDQGC
 WTYGNIRVSVEDVTVLVDYTVRKFCIQQVGDMTNRKPQRLITQFHFTSWPDFGVPFTPIGMLKFLKKVKA
 CNPQYAGAIVVHCSAGVGRGTGFVIDAMLDMHHTERKVDVYGFVSRIARQCQMVQTDQMYYFIIYQALL
 EHYLYGDTELEVTSLETHLQKIYNKIPGTSNNGLEEEFFKLLTSIKIQNDKMRTGNL PANMKKNRVLQIIP
 YEFNRVIVKVRGEENTDYVNASFIDGYRQKDSYIASQGPLLHTIEDFWRMIWEWKSCSIVMLTELEERG
 QEKCAQYWPSDGLVSYGDITVELKKEEESYTVRDLLVTNTRENKSRQIRQFHFGWPEVGIPSDGKGM
 ISIIAAVQKQQQSGNHPITVHCSAGAGRTGTFALSTVLERVKAEGILDVFQTVKSLRLQRPHMVQGLE
 QYEFYKVVQEYIDAFSDYANFK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	87.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_543031
Locus ID:	5786
UniProt ID:	P18433
RefSeq Size:	3153
Cytogenetics:	20p13
RefSeq ORF:	2379
Synonyms:	HEPTP; HLPR; HPTPA; HPTPalpha; LRP; PTPA; PTPRL2; R-PTP-alpha; RPTPA
Summary:	The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP contains an extracellular domain, a single transmembrane segment and two tandem intracytoplasmic catalytic domains, and thus represents a receptor-type PTP. This PTP has been shown to dephosphorylate and activate Src family tyrosine kinases, and is implicated in the regulation of integrin signaling, cell adhesion and proliferation. Three alternatively spliced variants of this gene, which encode two distinct isoforms, have been reported. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Phosphatase, Transmembrane

Product images:



Coomassie blue staining of purified PTPRA protein (Cat# [TP304392]). The protein was produced from HEK293T cells transfected with PTPRA cDNA clone (Cat# [RC204392]) using MegaTran 2.0 (Cat# [TT210002]).