

Product datasheet for **TP304422L**


SMR3B (NM_006685) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens submaxillary gland androgen regulated protein 3B (SMR3B), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204422 protein sequence Red =Cloning site Green =Tags(s) MKSLTWILGLWALAACFTPGESQGRGPYPYPPGPLAPPQPFPGFVPPPPPPYPGGRIPPPPPAPYGGP IFPPPPQP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_006676</u>
Locus ID:	10879
UniProt ID:	<u>P02814</u>
RefSeq Size:	798


[View online »](#)

Cytogenetics: 4q13.3
RefSeq ORF: 237
Synonyms: P-B; PBII; PRL3; PROL3; SMR1B
Protein Families: Secreted Protein

Product images:


Coomassie blue staining of purified SMR3B protein (Cat# [TP304422]). The protein was produced from HEK293T cells transfected with SMR3B cDNA clone (Cat# [RC204422]) using MegaTran 2.0 (Cat# [TT210002]).