

Product datasheet for TP304459

GAGE7 (NM_021123) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human G antigen 7 (GAGE7), full length, with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204459 protein sequence Red =Cloning site Green =Tags(s)
	MSWRGRSTYYWPRPRRYVQPPEVIGPMRPEQFSDEVEPATPEEGEPATQRQDPAAAQEGEDEGASAGQGP KPEADSQEQGHPQTGCECEDGPDGQEMDPPNPPEVKTPPEEGEKQSQC
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	Myc-DDK
Predicted MW:	12.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_066946
Locus ID:	2579
UniProt ID:	Q76087 , P0CL80 , P0CL81 , P0CL82
RefSeq Size:	498
Cytogenetics:	Xp11.4-p11.2
RefSeq ORF:	351



[View online »](#)

Synonyms: AL4; CT4.7; GAGE-7

Product images:



Coomassie blue staining of purified GAGE7 protein (Cat# TP304459). The protein was produced from HEK293T cells transfected with GAGE7 cDNA clone (Cat# [RC204459]) using MegaTran 2.0 (Cat# [TT210002]).