

Product datasheet for TP304515

ATCAY (NM_033064) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ataxia, cerebellar, Cayman type (ATCAY), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204515 protein sequence Red =Cloning site Green =Tags(s)
	<p>MGTTEATLRMENV DVKKEWQDEDLPRPLPEETGV ELLGSPVEDTSSPPNTLNFNGAHRKRKTLVAPEINI SLDQSEGSLLSDDFLDTPDDLINVDIETPDETSLEFLGNGNELEWEDDTPVATAKNMPGDSADLFGD GTTEDGSAANGRLWRTVIIGEQEHRIDLHMIRPYMKVWTHGGYYGEGLNAIIVFAACFLPDSSLPDYHYI MENFLYVISSLELLVAEDYMIVYLN GATPRRRMPGIGWLKCKYQMIDRRLRKLNKSLIIVHPSWFIRTV LAISRPFISVKFINKIQYVHSLEDLEQLIPMEHVQIPDCVLQYEEERL KARRESARPQPEFVLP RSEEK P EVAPVENRSALVSEDQETSMS</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	41.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_149053
Locus ID:	85300



[View online »](#)

UniProt ID:	Q86WG3 , A0A0S2Z5T8
RefSeq Size:	5049
Cytogenetics:	19p13.3
RefSeq ORF:	1113
Synonyms:	BNIP-H; CLAC
Summary:	This gene encodes a neuron-restricted protein that contains a CRAL-TRIO motif common to proteins that bind small lipophilic molecules. Mutations in this gene are associated with cerebellar ataxia, Cayman type. [provided by RefSeq, Jul 2008]
Protein Families:	Transmembrane

Product images:



Coomassie blue staining of purified ATCAY protein (Cat# TP304515). The protein was produced from HEK293T cells transfected with ATCAY cDNA clone (Cat# [RC204515]) using MegaTran 2.0 (Cat# [TT210002]).