

Product datasheet for TP304524M

CTRP6 (C1QTNF6) (NM_031910) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human C1q and tumor necrosis factor related protein 6 (C1QTNF6), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204524 protein sequence Red =Cloning site Green =Tags(s)
	<p>MQWLRVRESPGEATGHRVTMVTAALGPVWAALLLFLLMCEIPMVELTFDRAVASGCQRCCDSEDPLDPAH VSSASSGRPHALPEIRPYINITILKGDKGDPGPMGLPGYMGREGPQGEPGPQGSKGDKGEMGSPGAPCQ KRFFAFSVGRKTALHSGEDFQTLLEFRVFNLDGCFDMATGQFAAPLRGIYFFSLNVHSWNYKETYVHIM HNQKEAVILYAQPERSIMQSQSVMLDLAYGDRVWVRLFKRQRENAIYSNDFDTYITFSGHLIKAEDD</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	30.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_114116</u>
Locus ID:	114904



[View online »](#)

UniProt ID:	Q9BXI9 , A0A024R1J0
RefSeq Size:	2954
Cytogenetics:	22q12.3
RefSeq ORF:	834
Synonyms:	CTFP6; CTRP6; ZACRP6
Protein Families:	Secreted Protein, Transmembrane

Product images:



Coomassie blue staining of purified C1QTNF6 protein (Cat# [TP304524]). The protein was produced from HEK293T cells transfected with C1QTNF6 cDNA clone (Cat# [RC204524]) using MegaTran 2.0 (Cat# [TT210002]).