

## Product datasheet for **TP304538L**

### FGF21 (NM\_019113) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human fibroblast growth factor 21 (FGF21), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204538 representing NM_019113 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MDSDETGFHESGLWWSVLGALLLGGACQAHPIPDSSPLLQFGGQVRQRYLYTDDAQQTEAHLEIREDGTVG GAADQSPESLLQLKALKPGVIQILGVKTSRFLCQRPDGALYGSLHFDPEACSFRELLLEDGYNVYQSEAH GLPLHLPGNKSPHRDPAPRGPAPRFLPLPGLPPAPPEPPGILAPQPPDVGSSDPLSMVGPSQGRSPSYAS  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	19.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_061986</a>
Locus ID:	26291
UniProt ID:	<a href="#">Q9NSA1</a>
RefSeq Size:	940



[View online »](#)

Cytogenetics: 19q13.33


RefSeq ORF: 627

**Summary:** This gene encodes a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities and are involved in a variety of biological processes. This protein is a secreted endocrine factor that functions as a major metabolic regulator. The encoded protein stimulates the uptake of glucose in adipose tissue. [provided by RefSeq, Mar 2016]

**Protein Families:** Adult stem cells, Embryonic stem cells, ES Cell Differentiation/IPS, Secreted Protein

**Protein Pathways:** MAPK signaling pathway, Melanoma, Pathways in cancer, Regulation of actin cytoskeleton

### Product images:



Coomassie blue staining of purified FGF21 protein (Cat# [TP304538]). The protein was produced from HEK293T cells transfected with FGF21 cDNA clone (Cat# [RC204538]) using MegaTran 2.0 (Cat# [TT210002]).