

Product datasheet for **TP304560**

C11orf68 (NM_031450) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 11 open reading frame 68 (C11orf68), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204560 protein sequence Red =Cloning site Green =Tags(s)

MEPGEELEEEGSPGGREDGFTAETHLAAEAMAADMDPWLVFDARTTPATELDAWLAKYPPSQVTRYGDPGS
PNSEPVGWIAVYGQGYSPNSGDVQGLQAWEALQTSGRPITPGTLRQLAITHHVLVSGKWLMLHAPGFKLD
HAWAGIARAVVEGRLQVAKVSPRAKEGGRQVICVYTDFTDRLGVLEADSAIRAAGIKCLLTYKPDVYTY
LGIYRANRWHLCPITYESRFQLGGSARGSRVLDRAANNVELT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	27.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_113638.1
Locus ID:	83638



[View online »](#)

UniProt ID: [Q9H3H3](#), [Q9H3H3-2](#)
RefSeq Size: 1568
Cytogenetics: 11q13.1
RefSeq ORF: 753
Synonyms: BLES03; P5326

Product images:



Coomassie blue staining of purified C11orf68 protein (Cat# TP304560). The protein was produced from HEK293T cells transfected with C11orf68 cDNA clone (Cat# [RC204560]) using MegaTran 2.0 (Cat# [TT210002]).