

## Product datasheet for **TP304624**

### **GSTA3 (NM\_000847) Human Recombinant Protein**

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glutathione S-transferase alpha 3 (GSTA3), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204624 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MAGKPKLHYFNGRGRMEPIRWLLAAAGVEFEEKFIGSAEDLGKLRNDGSLMFQQVPMVEIDGIKLVQTRA ILNYIASKYNLYGKDIKERALIDMYTEGMADLNEMILLPLCRPEEKDAKIALIKEKTKSRYFPAFEKVL QSHGQDYLVGNKLSRADISLVLLYVEELDSSLISNFPLLKALKTRISNLPTVKKFLQPGSPRKPPADA KALEEARKIFRF
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	25.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_000838</a></u>
Locus ID:	2940
UniProt ID:	<u><a href="#">Q16772</a></u>



[View online »](#)

RefSeq Size:	915
Cytogenetics:	6p12.2
RefSeq ORF:	666
Synonyms:	GSTA3-3; GTA3
Summary:	<p>Cytosolic and membrane-bound forms of glutathione S-transferase are encoded by two distinct supergene families. These enzymes are involved in cellular defense against toxic, carcinogenic, and pharmacologically active electrophilic compounds. At present, eight distinct classes of the soluble cytoplasmic mammalian glutathione S-transferases have been identified: alpha, kappa, mu, omega, pi, sigma, theta and zeta. This gene encodes a glutathione S-transferase belonging to the alpha class genes that are located in a cluster mapped to chromosome 6. Genes of the alpha class are highly related and encode enzymes with glutathione peroxidase activity. However, during evolution, this alpha class gene diverged accumulating mutations in the active site that resulted in differences in substrate specificity and catalytic activity. The enzyme encoded by this gene catalyzes the double bond isomerization of precursors for progesterone and testosterone during the biosynthesis of steroid hormones. An additional transcript variant has been identified, but its full length sequence has not been determined. [provided by RefSeq, Jul 2008]</p>
Protein Pathways:	Drug metabolism - cytochrome P450, Glutathione metabolism, Metabolism of xenobiotics by cytochrome P450

### Product images:



Coomassie blue staining of purified GSTA3 protein (Cat# TP304624). The protein was produced from HEK293T cells transfected with GSTA3 cDNA clone (Cat# [RC204624]) using MegaTran 2.0 (Cat# [TT210002]).