

## Product datasheet for **TP304658L**

### CDKN2AIP (NM\_017632) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human CDKN2A interacting protein (CDKN2AIP), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204658 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MAQEVSEYLSQNPRVAAWVEALRCDGETDKHWRHRRDFLLRNAGDLAPAGGAASASTDEAADAESGTR NR</p> <p>QLQQLISFSMAWANHVFLGCRYPQKVMKILSMAEGIKVTDAPTYTTRDELVAKVKKRGISSSNEGVEEP SKKRVIIEGKNSSAVEQDHAKTSAKTERASAQQENSSTCIGSAIKSESGNSARSSGISSQNSSTSDGDRSV SSQSSSSVSSQVTTAGSGKASEAEAPDKHGSSFVSLKSSVNSHMTQSTDSRQQSGSPKKSALLEGSSASA SQSSSEIEVPLLGSSESEVELPLLSSKPSSETASSGLTSKTSSEASVSSSVAKNSSSSGTSLTPKSSS STNTSLLTSKSTSQVAASLLASKSSSQTSGSLVSKSTSLASVSQASKSSSQTSTSQLPSKSTSQSSESS VKFSCKLTNEDVKQKQPFFNRLYKTVAWKLVAVGGFSPNVNHGELLNAAIEALKATLDVFFVPLKELADL PQNKSSQESIVCELCKSVYLGTCGKSKENAKAVASREALKFLKKKVVVKICKRKYRGSEIEDLVLLD EESRPVNLPPALKHPQELL</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	60.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<u>NP_060102</u>
<b>Locus ID:</b>	55602
<b>UniProt ID:</b>	<u>Q9NXV6</u>
<b>RefSeq Size:</b>	2391
<b>Cytogenetics:</b>	4q35.1
<b>RefSeq ORF:</b>	1737
<b>Synonyms:</b>	CARF
<b>Summary:</b>	The protein encoded by this gene regulates the DNA damage response through several different signaling pathways. One such pathway is the p53-HDM2-p21(WAF1) pathway, which is critical to the DNA damage response. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2015]

### Product images:



Coomassie blue staining of purified CDKN2AIP protein (Cat# [TP304658]). The protein was produced from HEK293T cells transfected with CDKN2AIP cDNA clone (Cat# [RC204658]) using MegaTran 2.0 (Cat# [TT210002]).