

Product datasheet for **TP304697**

SLPI (NM_003064) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human secretory leukocyte peptidase inhibitor (SLPI), 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC204697 protein sequence
Red=Cloning site **Green**=Tags(s)

MKSSGLFPFLVLLALGTLAPWAVEGSGKSFKAGVCPPKKSQAQLRYKKPECQSDWQCPGKKRCCPDTCGI
KCLDPVDTPNPTRRRKPGKCPVTYGQCLMLNPPNFCEMDGQCKRDLKCCMGMCVKSCVSPVKA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 11.7 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_003055](#)

Locus ID: 6590

UniProt ID: [P03973](#)

RefSeq Size: 625

Cytogenetics: 20q13.12



[View online »](#)


RefSeq ORF: 396

Synonyms: ALK1; ALP; BLPI; HUSI; HUSI-I; MPI; WAP4; WFDC4

Summary: This gene encodes a secreted inhibitor which protects epithelial tissues from serine proteases. It is found in various secretions including seminal plasma, cervical mucus, and bronchial secretions, and has affinity for trypsin, leukocyte elastase, and cathepsin G. Its inhibitory effect contributes to the immune response by protecting epithelial surfaces from attack by endogenous proteolytic enzymes. This antimicrobial protein has antibacterial, antifungal and antiviral activity. [provided by RefSeq, Nov 2014]

Protein Families: Secreted Protein

Product images:



Coomassie blue staining of purified SLPI protein (Cat# TP304697). The protein was produced from HEK293T cells transfected with SLPI cDNA clone (Cat# [RC204697]) using MegaTran 2.0 (Cat# [TT210002]).