

Product datasheet for **TP304733M**

Thyroxine Binding Globulin (SERPINA7) (NM_000354) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 7 (SERPINA7), 100 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC204733 protein sequence
Red=Cloning site **Green**=Tags(s)

MSPFLYLVLVLLVGLHATHCASPEGKVTACHSSQP NATLYKMSSINADFAFNLYRRFTVETPDKNIFFSP
VSISAALVMLSFGACCSTQTEIVETLGFNLTDTPMVEIQHGFQHLICSLNFPKKELELQIGNALFIGKHL
KPLAKFLNDVKTYETEVFSTDFSNISA AKQEINSHVEMQTKGKVVGLIQDLKPNTIMVLVNYIHFKAQW
ANPFDPSKTEDSSSFLIDKTTTQVPMMHQMEQYYHLVDMELNCTVLQMDYSKNALALFVLPKEGQMESV
EAAMSSKTLKKWNRLQKGVVDFVFPKFSISATYDLGATLLKMGIQHAYSENADFSGLTEDNGLKLSNAA
HKAVLHIGEKGTEAAAVPEVELSDQPENTFLHPHIIQIDRSFMLLILERSTRSILFLGKVVNPT EA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 46.1 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_000345](#)




[View online »](#)

Locus ID: 6906
UniProt ID: [P05543](#)
RefSeq Size: 1600
Cytogenetics: Xq22.3
RefSeq ORF: 1245
Synonyms: TBG; TBGQTL

Summary: There are three proteins including thyroxine-binding globulin (TBG), transthyretin and albumin responsible for carrying the thyroid hormones thyroxine (T4) and 3,5,3'-triiodothyronine (T3) in the bloodstream. This gene encodes the major thyroid hormone transport protein, TBG, in serum. It belongs to the serpin family in genomics, but the protein has no inhibitory function like many other members of the serpin family. Mutations in this gene result in TGB deficiency, which has been classified as partial deficiency, complete deficiency, and excess, based on the level of serum TBG. Alternatively spliced transcript variants encoding different isoforms have been found, but the full-length nature of these variants has not been determined.[provided by RefSeq, Jun 2012]

Protein Families: Druggable Genome, Secreted Protein

Product images:



Coomassie blue staining of purified SERPINA7 protein (Cat# [TP304733]). The protein was produced from HEK293T cells transfected with SERPINA7 cDNA clone (Cat# [RC204733]) using MegaTran 2.0 (Cat# [TT210002]).