

Product datasheet for TP304760M

SHKBP1 (NM_138392) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human SH3KBP1 binding protein 1 (SHKBP1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204760 protein sequence Red =Cloning site Green =Tags(s)
	<p>MAAAATAAEGVPSRGPPEGEVIHLNVGGKRFSTSRQTLTWIPDSFFSLLSGRISTLKDETGAIFIDRDPT VFAPILNFLRTKELDPRGVHGSSLLHEAQFYGLTPLVRRQLREELDRSSCGNVLFNGYLPPPVPVKRR NRHSLVGPQQLGGRPAPVRRSNTMPPNLGNAGLLGRMLDEKTPPSPSGQPEEPMVRLVCGHHNWIAYAV TQFLVCYRLKEASGWQLVFSSPRLDWPIERLALTARVHGGALGEHDKMVAATGSEILLWALQAEGGGSE IGVFHLGVPVEALFFVGNQLIATSHTGRIGVWNAVTKHWQVQEVQPITSYDAAGSLLLGCNNGSIYYVD VQKFPLRMKDNDLLVSELYRDPADGVTALSVYLPKTSDSGNWIEIAYGTSSGGVRVIVQHPETVGGSP QLFQFTFVHRSPVTKIMLSEKHLISVCADNNHVRTWSVTRFRGMISTQPGSTPLASFKILALESADGHGG CSAGNDIGPYGERDDQLVFIQKVPSASQLFVRLSSTGQRVCSVRSVDGSPPTAFTVLECEGSRRLGSRP RRYLLTGQANGSLAMWDLTTAMDGLGQAPAGGLTEQELMEQLEHCELAPPAPSAPSWGCLPSPSPRISLT SLHSASSNTSLSGHRGSPSPQAARRRRGGGSFVERCQELVRSRGPDLRRPPTPAPWPSSGLGTPLTPPKM KLNETSF</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	76.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_612401</u>
Locus ID:	92799
UniProt ID:	<u>Q8TBC3</u>
RefSeq Size:	2379
Cytogenetics:	19q13.2
RefSeq ORF:	2121
Synonyms:	PP203; Sb1
Summary:	Inhibits CBL-SH3KBP1 complex mediated down-regulation of EGFR signaling by sequestration of SH3KBP1. Binds to SH3KBP1 and prevents its interaction with CBL and inhibits translocation of SH3KBP1 to EGFR containing vesicles upon EGF stimulation.[UniProtKB/Swiss-Prot Function]
Protein Families:	Ion Channels: Other

Product images:

Coomassie blue staining of purified SHKBP1 protein (Cat# [TP304760]). The protein was produced from HEK293T cells transfected with SHKBP1 cDNA clone (Cat# [RC204760]) using MegaTran 2.0 (Cat# [TT210002]).