

Product datasheet for **TP304762L**

RIL (PDLIM4) (NM_003687) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human PDZ and LIM domain 4 (PDLIM4), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204762 protein sequence Red =Cloning site Green =Tags(s)
	<p>MPHSVTLRGPSWPWFRLVGGGRDFSAPLTISR VHAGSKAALAALCPGDLIQAINGESTELMTHLEAQNRIK GCHDHLTLSVSRPEGRSWPSAPDDSKAQAHRIHIDPEIQDGSPTTSRGPSTGTGPEDGRPSLGSPYGQP PRFPVPHNGSSEATLPAQMSTLHVSPPSADPARGLPRSRDCRVDLGSEVYRMLREPAEPVAAEPKQSGS FRYLQGMLEAGEGGDWPGPGGPRNLKPTASKLGLAPLSGLQGLPECTRCGHGIVGTIVKARDKLYHPECFM CSDCGLNLKQRGYFFLDERLYCESHAKARVKPPEGYDVVAVYPNAKVELV</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	35.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_003678</u>
Locus ID:	8572



[View online »](#)

UniProt ID: [P50479](#)

RefSeq Size: 2307


Cytogenetics: 5q31.1

RefSeq ORF: 990

Synonyms: RIL

Summary: This gene encodes a protein which may be involved in bone development. Mutations in this gene are associated with susceptibility to osteoporosis. [provided by RefSeq, Nov 2009]

Product images:



Coomassie blue staining of purified PDLIM4 protein (Cat# [TP304762]). The protein was produced from HEK293T cells transfected with PDLIM4 cDNA clone (Cat# [RC204762]) using MegaTran 2.0 (Cat# [TT210002]).