

Product datasheet for TP304784L

SERF2 (NM_001018108) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human small EDRK-rich factor 2 (SERF2), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204784 protein sequence Red =Cloning site Green =Tags(s) MTRGNQRELARQKNMKQSDSVK GKRRDDGLSAAARKQRDSEIMQQKQKKANEKKEEPPK TR TRPLE QKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	6.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001018118</u>
Locus ID:	10169
UniProt ID:	<u>P84101</u>
RefSeq Size:	3048
Cytogenetics:	15q15.3


[View online »](#)

RefSeq ORF:	177
Synonyms:	4F5REL; FAM2C; H4F5REL; HsT17089
Summary:	Positive regulator of amyloid protein aggregation and proteotoxicity (PubMed:20723760). Induces conformational changes in amyloid proteins, such as HTT, driving them into compact formations preceding the formation of aggregates (PubMed:20723760).[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified SERF2 protein (Cat# [TP304784]). The protein was produced from HEK293T cells transfected with SERF2 cDNA clone (Cat# [RC204784]) using MegaTran 2.0 (Cat# [TT210002]).