

## Product datasheet for TP304786

### MRPL46 (NM\_022163) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human mitochondrial ribosomal protein L46 (MRPL46), nuclear gene encoding mitochondrial protein, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204786 protein sequence Red=Cloning site Green=Tags(s)

MAAPVRRTLLGVAGGWRRFERLWAGSLSSRSLALAAAPSSNGSPWRLLGALCLQRPPVWSKPLTPLQEEM  
ASLLQQIEIERSLYSDHELALDENQRLAKKKADLHDEEDEQDILLAQDLEDMWEQKFLQFKLGARITEA  
DEKNDRTSLNRKLDNRNLVLLVREKFGDQDVWILPQAEWQPGETLRGTAERTLATLSENMEAKFLGNAPC  
GHYTFKFPQAMRTESNLGAKVFFKALLLTGDFSQAGNKGHHVWVKDELGDYLPKYLAQVRRFVSDL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	31.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_071446</a></u>
Locus ID:	26589



[View online »](#)

UniProt ID: [Q9H2W6](#)

RefSeq Size: 1020

Cytogenetics: 15q25.3

RefSeq ORF: 837

Synonyms: C15orf4; LIECG2; P2ECSL

**Summary:** Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. [provided by RefSeq, Jul 2008]

### Product images:

□

Coomassie blue staining of purified MRPL46 protein (Cat# TP304786). The protein was produced from HEK293T cells transfected with MRPL46 cDNA clone (Cat# [RC204786]) using MegaTran 2.0 (Cat# [TT210002]).