

## Product datasheet for TP304793L

### PLAC1 (NM\_021796) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human placenta-specific 1 (PLAC1), with C-terminal MYC/DDK tag, expressed in HEK293 cells, 1 mg
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	>RC204793 representing NM_021796 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MKVFKFIGLMILLTSAFSAGSGQSPMTVLC SIDWFMVTVHPFMLNNDVCVHFHELHLGLGCPPNHVQPHA YQFTYRVTECGIRAKAVSQDMVIYSTEIHYSKGTPSKVFVIVS CAAPQKSPWLTKPCSMRVASKRATA QKDEKCYEVFSLSQSSQRPNDCPPCVFSEEHTQVPCHQAGAQAQPLQPSHFLDISEDWSLHTDDMIG SM
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	23.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_068568</a></u>
Locus ID:	10761



[View online »](#)

UniProt ID:	<a href="#">Q9HBJ0</a>
RefSeq Size:	1126
Cytogenetics:	Xq26.3
RefSeq ORF:	636
Synonyms:	CT92; OOSP2B; OOSP2L
Summary:	May play a role in placental development.[UniProtKB/Swiss-Prot Function]
Protein Families:	Secreted Protein

### Product images:



Coomassie blue staining of purified PLAC1 protein (Cat# [TP304793]). The protein was produced from HEK293T cells transfected with PLAC1 cDNA clone (Cat# [RC204793]) using MegaTran 2.0 (Cat# [TT210002]).