

# **Product datasheet for TP304804**

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

#### BAP31 (BCAP31) (NM 005745) Human Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human B-cell receptor-associated protein 31 (BCAP31), transcript

variant 2, 20 µg

Species: Human
Expression Host: HEK293T

**Expression cDNA Clone** >RC204804 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s)

MSLQWTAVATFLYAEVFVVLLLCIPFISPKRWQKIFKSRLVELLVSYGNTFFVVLIVILVLLVIDAVREI REYDDVTEKVNLQNNPGAMEHFHMKLFRAQRNLYIAGFSLLLSFLLRRLVTLISQQATLLASNEAFKKQA ESASEAAKKYMEENDQLKKGAAVDGGKLDVGNAEVKLEEENRSLKADLQKLKDELASTKQKLEKAENQVL

AMRKQSEGLTKEYDRLLEEHAKLQAAVDGPMDKKEE

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV** 

Tag: C-Myc/DDK

**Predicted MW:** 27.8 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** NP 005736

**Locus ID:** 10134



#### BAP31 (BCAP31) (NM\_005745) Human Recombinant Protein - TP304804

UniProt ID: P51572
RefSeq Size: 1417
Cytogenetics: Xq28
RefSeq ORF: 738

Synonyms: 6C6-AG; BAP31; CDM; DDCH; DXS1357E


Summary: This gene encodes a member of the B-cell receptor associated protein 31 superfamily. The

encoded protein is a multi-pass transmembrane protein of the endoplasmic reticulum that is involved in the anterograde transport of membrane proteins from the endoplasmic reticulum to the Golgi and in caspase 8-mediated apoptosis. Microdeletions in this gene are associated

with contiguous ABCD1/DXS1375E deletion syndrome (CADDS), a neonatal disorder. Alternative splicing of this gene results in multiple transcript variants. Two related pseudogenes have been identified on chromosome 16. [provided by RefSeq, Jan 2012]

**Protein Families:** Druggable Genome, Transmembrane

## **Product images:**



Coomassie blue staining of purified BCAP31 protein (Cat# TP304804). The protein was produced from HEK293T cells transfected with BCAP31 cDNA clone (Cat# [RC204804]) using

MegaTran 2.0 (Cat# [TT210002]).