

Product datasheet for TP304810M

CA12 (NM_001218) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human carbonic anhydrase XII (CA12), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204810 protein sequence Red =Cloning site Green =Tags(s)

MPRRSLHAAAVLLLVLKEQPSSPAPVNGSKWTFYFGPDGENSWSKYPSCGGLLQSPIDLHSDILQYDAS
LTPLEFQGYNLSANKQFLLTNNHGSVKLNLPSPDMHIQGLQSRYSATQLHLHWGNPNDPHGSEHTVSGQHF
AAELHIVHYNSDLYPDASTASNKSEGLAVLAVLIEMGSFNPSYDKIFSHLQHVKYKGQEA FVPGFNIEEL
LPERTA EYRYRGSLLTPPCNPTVLWTVFRNPVQISQEQLLALETALYCTHMDDPSPREMINNFRQVQKF
DERLVYTSFSQVQVCTAAGLSLGIILSLALAGILGICIVVVSIWLFRRRSIKKGDNKGVIYKPKMET
EAHA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	36.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001209
Locus ID:	771



[View online »](#)

UniProt ID: [Q43570](#)

RefSeq Size: 4209

Cytogenetics: 15q22.2

RefSeq ORF: 1062


Synonyms: CA-XII; CAXII; HsT18816; T18816

Summary: Carbonic anhydrases (CAs) are a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. They participate in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva, and gastric acid. This gene product is a type I membrane protein that is highly expressed in normal tissues, such as kidney, colon and pancreas, and has been found to be overexpressed in 10% of clear cell renal carcinomas. Three transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jun 2014]

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Nitrogen metabolism

Product images:



Coomassie blue staining of purified CA12 protein (Cat# [TP304810]). The protein was produced from HEK293T cells transfected with CA12 cDNA clone (Cat# [RC204810]) using MegaTran 2.0 (Cat# [TT210002]).