

## Product datasheet for **TP304819**

### **RPB11 (POLR2J) (NM\_006234) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human polymerase (RNA) II (DNA directed) polypeptide J, 13.3kDa (POLR2J), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204819 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MNAPPAFESFLLFEGEKKITINKDTPVNAFLFTINKEDHTLGNIIKSQLLKDPQVLFAGYKVPHPLEHK  
IIRVQTTPDYSPQEAFTNAITDLISELSLLEERFRVAIKDKQEGIE

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**


Tag:	C-Myc/DDK
Predicted MW:	13.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_006225</a>
Locus ID:	5439
UniProt ID:	<a href="#">P52435</a>
RefSeq Size:	991



[View online »](#)

Cytogenetics:	7q22.1
RefSeq ORF:	351
Synonyms:	hRPB14; POLR2J1; RPB11; RPB11A; RPB11m
Summary:	This gene encodes a subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. The product of this gene exists as a heterodimer with another polymerase subunit; together they form a core subassembly unit of the polymerase. Two similar genes are located nearby on chromosome 7q22.1 and a pseudogene is found on chromosome 7p13. [provided by RefSeq, Jul 2008]
Protein Families:	Transcription Factors
Protein Pathways:	Huntington's disease, Metabolic pathways, Purine metabolism, Pyrimidine metabolism, RNA polymerase

### Product images:



Coomassie blue staining of purified POLR2J protein (Cat# TP304819). The protein was produced from HEK293T cells transfected with POLR2J cDNA clone (Cat# [RC204819]) using MegaTran 2.0 (Cat# [TT210002]).