

Product datasheet for **TP304823L**

MIRO2 (RHOT2) (NM_138769) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ras homolog gene family, member T2 (RHOT2), nuclear gene encoding mitochondrial protein, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204823 protein sequence Red =Cloning site Green =Tags(s)

MRRDVRILLLGEAQVGKTSLILSLVGEEFPPEEVPRAEEITIPADVTPKVPHTHIVDYSEAEQTDEELRE
EIHKANVVCVYDVSEETIEKIRTKWIPLVNGGTTQGPRVPIILVGNKSDLRSGSSMEAVLPIMSQFPE
IETCVECSAKNLRNISELFYQAQKAVLHPTAPLYDPEAKQLRPACAQALTRIFRLSDQDLQALSDEELN
AFQKSCFGHPLAPQALEDVKTVCVCRNVAGGVREDRLTLDGFLFLNTLFIQRGRHETTWTILRRFGYS DAL
ELTADYLSPLIHVPPGCSTELNHLGYQFVQRVFEKHDQDRDGALSPVELQSLFSVFPAAWGPPELPRTVR
TEAGRLPLHGYLCQWTLVTYLDVRSCLGHLGYLGYPTLCEQDQAHAITVTRKRLDQEKGQTQRSVLLCK
VVGARGVGKSAFLQAFGLRGLGHQDTREQPPGYAIDTVQVNGQEKYLILCEVGTDGLLATS DATCDVAC
LMFDGSDPKSFAHCASVYKHHYMDGQTPCLFVSSKADLPEGVAVSGPSPAEFCRKHRLPAPVPFSCAGPA
EPSTTIFTQLATMAAFPHLVHAELHPSSFVLRGLLGVGAAVA AVLSFSLYRVLVKSQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	67.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_620124
Locus ID:	89941
UniProt ID:	Q8IXI1
RefSeq Size:	2556
Cytogenetics:	16p13.3
RefSeq ORF:	1854
Synonyms:	ARHT2; C16orf39; MIRO-2; MIRO2; RASL
Summary:	This gene encodes a member of the Rho family of GTPases. The encoded protein is localized to the outer mitochondrial membrane and plays a role in mitochondrial trafficking and fusion-fission dynamics. [provided by RefSeq, Nov 2011]
Protein Families:	Transmembrane

Product images:

Coomassie blue staining of purified RHOT2 protein (Cat# [TP304823]). The protein was produced from HEK293T cells transfected with RHOT2 cDNA clone (Cat# [RC204823]) using MegaTran 2.0 (Cat# [TT210002]).