

## Product datasheet for TP304825M

### NMNAT1 (NM\_022787) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nicotinamide nucleotide adenylyltransferase 1 (NMNAT1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204825 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MENSEKTEVLLACGSFNPITNMHLRLFELAKDYMNGTGRYTVWKGIIISPVGDAYKKKGLIPAYHRVIMA  
ELATKNSKWVEVDTWESLQKEWKETLKVLRHHQEKLASDCDHQQNSPTLERPGRKRKWTETQDSSQKKS  
LEPKTKAVPKVKLLCGADLLESFAVPLWKSEDITQIVANYGLICVTRAGNDAQKFIYESDVLWKHRSNI  
HVVNEWIANDISSTKIRRALRRGQSIRYLVPDLVQEYIEKHNLYSSESEDRNAGVILAPLQRNTAEAKT

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	31.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_073624</a>
Locus ID:	64802



[View online »](#)

UniProt ID: [Q9HAN9](#), [A0A024R4E1](#)

RefSeq Size: 3781

Cytogenetics: 1p36.22


RefSeq ORF: 837

Synonyms: LCA9; NMNAT; PNAT1; SHILCA

**Summary:** This gene encodes an enzyme which catalyzes a key step in the biosynthesis of nicotinamide adenine dinucleotide (NAD). The encoded enzyme is one of several nicotinamide nucleotide adenylyltransferases, and is specifically localized to the cell nucleus. Activity of this protein leads to the activation of a nuclear deacetylase that functions in the protection of damaged neurons. Mutations in this gene have been associated with Leber congenital amaurosis 9. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are located on chromosomes 1, 3, 4, 14, and 15. [provided by RefSeq, Jul 2014]

**Protein Pathways:** Metabolic pathways, Nicotinate and nicotinamide metabolism

### Product images:



Coomassie blue staining of purified NMNAT1 protein (Cat# [TP304825]). The protein was produced from HEK293T cells transfected with NMNAT1 cDNA clone (Cat# [RC204825]) using MegaTran 2.0 (Cat# [TT210002]).