

Product datasheet for **TP304828L**

YARS2 (NM_001040436) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tyrosyl-tRNA synthetase 2, mitochondrial (YARS2), nuclear gene encoding mitochondrial protein, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204828 protein sequence Red =Cloning site Green =Tags(s)

MAAPILRSFSWGRWSGTLNLSVLLPLGLRKAHSGAQGLLAAQKARGLFKDFFPETGTKIELPELFDRTGTA
SFPQTIYCGFDPTADSLHVGHLLALLGLFHLQRAGHNVIALVGGATARLGDPGRTKEREALETERVRAN
ARALRLGLEALAAHQQLFTDGRSWSFTVLDNSAWYQKQHLVDFLAAVGGHFRMGTLRSQSVQLRLKS
PEGMSLAEFFYQVLQAYDFYLLFQRYGCRVQLGGSDQLGNIMSGYEFINKLTGEDVFGITVPLITSTTGA
KLGKSAGNAVWLNDRKTSPELYQFFVRQPDDSVERYLKLFTFLPLPEIDHIMQLHVKEPERRGPQKRLA
AEVTKLVHGREGLDSAKRCTQALYHSSIDALEVMSDQELKELFKEAPFSEFFLDPGTSVLDTCRKANAIIP
DGPRGYRMITEGGVSINHQQVTNPESVLIVGQHILKNGLSLLKIGKRNFYIHKWLQL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	53 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq: [NP_001035526](#)

Locus ID: 51067

UniProt ID: [Q9Y2Z4](#)

RefSeq Size: 2173

Cytogenetics: 12p11.21

RefSeq ORF: 1431


Synonyms: CGI-04; MLASA2; MT-TYRRS; TYRRS

Summary: This gene encodes a mitochondrial protein that catalyzes the attachment of tyrosine to tRNA(Tyr). Mutations in this gene are associated with myopathy with lactic acidosis and sideroblastic anemia type 2 (MLASA2). [provided by RefSeq, Jan 2011]

Protein Families: Druggable Genome

Protein Pathways: Aminoacyl-tRNA biosynthesis

Product images:



Coomassie blue staining of purified YARS2 protein (Cat# [TP304828]). The protein was produced from HEK293T cells transfected with YARS2 cDNA clone (Cat# [RC204828]) using MegaTran 2.0 (Cat# [TT210002]).