

Product datasheet for **TP304856L**

CMPK1 (NM_016308) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cytidine monophosphate (UMP-CMP) kinase 1, cytosolic (CMPK1), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204856 protein sequence Red =Cloning site Green =Tags(s) MLSRCRSRLHVLGLSFLQTRRPILLCSPLMKPLVWFVLGGPGAGKGTQCARIVEKYGYTHLSAGELL RDERKNPDSQYGELIEKYIKEGKIVPVEITISLLKREMDQTMAANAQKNKFLIDGFPRNQDNLQGWNKTM DGKADVSFVLFDCNNEICIERCLERKSSGRSDDNRESLEKRIQTYLQSTKPIIDL YEEMGKVKKIDAS KSVDEVFDEVVQIFDKEG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	25.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	Enzyme activity (PMID: 26167664)
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057392



[View online »](#)


Locus ID:	51727
UniProt ID:	P30085
RefSeq Size:	2956
Cytogenetics:	1p33
RefSeq ORF:	684
Synonyms:	CK; CMK; CMPK; UMK; UMP-CMPK; UMPK

Summary: This gene encodes one of the enzymes required for cellular nucleic acid biosynthesis. This enzyme catalyzes the transfer of a phosphate group from ATP to CMP, UMP, or dCMP, to form the corresponding diphosphate nucleotide. Alternate splicing results in both coding and non-coding transcript variants. [provided by RefSeq, Feb 2012]

Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Pyrimidine metabolism

Product images:



Coomassie blue staining of purified CMPK1 protein (Cat# [TP304856]). The protein was produced from HEK293T cells transfected with CMPK1 cDNA clone (Cat# [RC204856]) using MegaTran 2.0 (Cat# [TT210002]).