

## Product datasheet for TP304872

### XPA (NM\_000380) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Human xeroderma pigmentosum, complementation group A (XPA), transcript variant 1, full length, with C-terminal MYC/DDK tag, expressed in HEK293 cells, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204872 protein sequence Red=Cloning site Green=Tags(s)
	<p>MAAADGALPEAAALEQPAELPASVRASIERKRQRALMLRQARLAARPYSTATAAAATGGMANVKAAPKIID TGGGFILEEEEEEQKIGKVWHQPGPVMEDYVICEECGKEFMDSYLMNHFDLPTCDNCRDADDKHKLIT KTEAKQEYLLKDCDLEKREPPLKFIVKKNPHHSQWGMKLYLKLQIVKRSLEVWGSQEALEEAKEVRQEN REKMKQKKFKKVKELRRAVRSSVWKRETIVHQHEYGPEENLEDDMYRKTCTMCGHELTYEKM</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	31.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_000371</a>
Locus ID:	7507



[View online »](#)

UniProt ID: [P23025](#)

RefSeq Size: 1491

Cytogenetics: 9q22.33

RefSeq ORF: 819


Synonyms: XP1; XPAC

**Summary:** This gene encodes a zinc finger protein plays a central role in nucleotide excision repair (NER), a specialized type of DNA repair. NER is responsible for repair of UV radiation-induced photoproducts and DNA adducts induced by chemical carcinogens and chemotherapeutic drugs. The encoded protein interacts with DNA and several NER proteins, acting as a scaffold to assemble the NER incision complex at sites of DNA damage. Mutations in this gene cause Xeroderma pigmentosum complementation group A (XP-A), an autosomal recessive skin disorder featuring hypersensitivity to sunlight and increased risk for skin cancer. [provided by RefSeq, Aug 2017]

**Protein Families:** Druggable Genome

**Protein Pathways:** Nucleotide excision repair

### Product images:



Coomassie blue staining of purified XPA protein (Cat# TP304872). The protein was produced from HEK293T cells transfected with XPA cDNA clone (Cat# [RC204872]) using MegaTran 2.0 (Cat# [TT210002]).