

Product datasheet for **TP304914**

FAM136A (NM_032822) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 136, member A (FAM136A), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204914 protein sequence Red =Cloning site Green =Tags(s)
	MAELQQLRVQEAVESMVKSLERENIRKMQGLMFRCSASCCEDSQASMKQVHQCIERCHVPLAQAQALV TS ELEKFQDRLARCTMHCNDKAKDSIDAGSKELQVKQQLDSCVTKCVDHMHLIPTMTKKMKEALLSIGK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	15.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_116211</u>
Locus ID:	84908
UniProt ID:	<u>Q96C01</u>



[View online »](#)


RefSeq Size: 1824

Cytogenetics: 2p13.3

RefSeq ORF: 414

Summary: This gene encodes a mitochondrially localized protein that is highly conserved across species. The gene is expressed in a variety of tissues including human lymphoblast cells and rat neurosensory epithelium of the crista ampullaris. A mutation in this gene has been associated with familial Meniere's disease, a chronic disorder of the inner ear. Several pseudogenes of this gene are found on other chromosomes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2016]

Product images:



Coomassie blue staining of purified FAM136A protein (Cat# TP304914). The protein was produced from HEK293T cells transfected with FAM136A cDNA clone (Cat# [RC204914]) using MegaTran 2.0 (Cat# [TT210002]).