

Product datasheet for **TP304931L**

LENG1 (NM_024316) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human leukocyte receptor cluster (LRC) member 1 (LENG1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204931 protein sequence Red =Cloning site Green =Tags(s)

MNILPKKSWHVRNKDNVARVRRDEAQAREEEKERERRVLLAQQEARTFLRKKARHQNSLPELEAAEAGA
PGSGPVDLFRLEEGKGVIRGNKEYEEERQEKERQEKALGILTYLGQSAEAQTQPPWYQLPPGRGGP
PPGPAPDEKIKSRLDPLREMQKHLGKKRQHGGDEGSRSRKEKEGSEKQRPKEPPSLDQLRAERLRREAAE
RSRAEALLARVQGRALQEGQPEDETDDRRRRYNSQFNPQLARRPRQQDPHLTH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	30.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_077292
Locus ID:	79165
UniProt ID:	Q96BZ8




[View online >](#)

RefSeq Size: 941

Cytogenetics: 19q13.42

RefSeq ORF: 792

Product images:



Coomassie blue staining of purified LENG1 protein (Cat# [TP304931]). The protein was produced from HEK293T cells transfected with LENG1 cDNA clone (Cat# [RC204931]) using MegaTran 2.0 (Cat# [TT210002]).