

Product datasheet for **TP304934L**

AMID (AIFM2) (NM_032797) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human apoptosis-inducing factor, mitochondrion-associated, 2 (AIFM2), nuclear gene encoding mitochondrial protein, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204934 protein sequence Red =Cloning site Green =Tags(s)

MGSQVSVESGALHWVIVGGGFGGIAAASQLQALNVPFMLVDMKDSFHNNVAALRASVETGFAKKTFSYS
VTFKDNFRQGLVVGIDLKNQMVLLQGGEALPFSHLILATGSTGPFPGKFNEVSSQQAIIQAYEDMVRQVQ
RSRFIVVGGGSAGVEMAAEIKTEYPEKEVTLIHSQVALADKELLPSVRQEVEKILLRKGVQLLSERVS
NLEELPLNEYREYIKVQTDKGTAVTNLVILCTGIKINSSAYRKAFESRLASSGALRVNEHLQVEGHSNV
YAIGDCADV RTPKMAYLAGLHANI AVANIVNSVKQRPLQAYKPGALTFLLSMGRNDG VGVQISGFYVGRML
VRLTKSRDLFVSTSWKTMRQSP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	40.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_116186</u>



[View online »](#)

Locus ID:	84883
UniProt ID:	Q9BRQ8
RefSeq Size:	3253
Cytogenetics:	10q22.1
RefSeq ORF:	1119
Synonyms:	AMID; FSP1; PRG3
Summary:	This gene encodes a flavoprotein oxidoreductase that binds single stranded DNA and is thought to contribute to apoptosis in the presence of bacterial and viral DNA. The expression of this gene is also found to be induced by tumor suppressor protein p53 in colon cancer cells. [provided by RefSeq, Nov 2010]
Protein Families:	Druggable Genome, Transmembrane

Product images:



Coomassie blue staining of purified AIFM2 protein (Cat# [TP304934]). The protein was produced from HEK293T cells transfected with AIFM2 cDNA clone (Cat# [RC204934]) using MegaTran 2.0 (Cat# [TT210002]).