

Product datasheet for **TP304945L**

FAM219B (NM_020447) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 15 open reading frame 17 (C15orf17), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC204945 protein sequence Red =Cloning site Green =Tags(s)
	MATAEPSGRALRLSTPGPRPSGARDRAPGAAGPPSGQIGNRALRLGERTPAAVEKRGPYMVTRAPSIQAK LQKHRDLAKAVLRRKGMGLGASPNRPDSSGKRsvkfnkgytalsqspdenlvslsdSDGELGSRySSGYS SAEQVNQDVSRQLLQDGYHLDEIPDDEDLDLIPPKPMASSTCSCCWCLGDSSSCTLQ
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	20.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_065180</u>
Locus ID:	57184
UniProt ID:	<u>Q5XKK7</u>
RefSeq Size:	3563



[View online »](#)

Cytogenetics: 15q24.1-q24.2
RefSeq ORF: 594
Synonyms: C15orf17
Protein Families: Druggable Genome

Product images:

Coomassie blue staining of purified FAM219B protein (Cat# [TP304945]). The protein was produced from HEK293T cells transfected with FAM219B cDNA clone (Cat# [RC204945]) using MegaTran 2.0 (Cat# [TT210002]).